

EVLT FIBERS BY NEULIGHT®

SINGLE-USE AND REUSABLE



NeuLight is pleased to announce the approval of our 10x reusable and single-use 600um fiber for use in the laser based varicose vein therapy market for North and South America.

Max. input power	100 W
Length	3 m
Numerical Aperture	0,37
Core Size	600 um
Coating Diameter	630 um
Buffer Diameter	750 um
Connector Type	SMA905
Wavelength	500 nm – 2200 nm
Usability	Sterile, Ready to use with open-port lasers
Product Codes	Single-Use: WF IR600HCT-3/SM-F 2M20Z1C Reusable: WF IR600HCT-3/SM-F 3N20Z3C

EVLT or Endovenous Laser Therapy is the minimally invasive treatment method for varicose veins through the use of laser energy. Varicose veins are swollen, twisted veins that develop in legs and can cause discomfort. During this procedure, the heat produced by the laser closes off the varicose veins from blood flow. This forces the blood to find healthier veins, causing varicose veins and their associated symptoms to disappear.

- NeuLight offers the highest quality fibers at reasonable prices.
- NeuLight EVLT fibers are manufactured by WEINERT Fiber Optics, an ISO13485 certified group with a global presence.
- High quality laser fibers compatible with standard SMA-905 fiber connectors on EVLT laser systems.
- Designed with bare tips and applicable with wavelengths ranging from 500 nm – 2200 nm delivering up to 100 Watts of operating power.
- NeuLight EVLT fibers are available from our trusted distribution partners in the United States and South America.
- Probes are sold sterile in boxes of 5 units.
- Single-Use and 10x Reusable versions available.
- Flat-tip design
- EZ deliver and protective sterile packaging system for ease-of-use and safety in the operating room.
- NeuLight EVLT Fibers have been used in countless minimally invasive EVLT procedures world-wide.
- Private labeling is available, Integration of RFID technology possible.



Subject to technical changes without notice.

Weinert Fiber Optics Inc. · 209 Bulifants Blvd. · Williamsburg, VA 23188, USA · +1 (757) 258-4805 · sales-americas@weinert-industries.com